Superlite Power Pen/ Color Pen

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Superlite Power Pen/ Color Pen

Instructions for use

Thank you for trusting our products!

Carefully read this Operating instructions manual before installing and operating the unit to use and maintain it successfully!

Save this Operations guide for your reference in the future.

CONTENTS

I. UNIT DESCRIPTION AND FUNCTION

II. SYMBOLS

III. SAFETY PRECAUTIONS

IV. TECHNICAL DATA

V. COMPLETE SET

VI. PREPARATION AND SEQUENCE OF OPERATION

VII. ROUTINE CARE AND MAINTENANCE

VIII. WARRANTY CONDITIONS

IX. SERVICE DATA

X. FAQ

XI. DECLARATION OF CONFORMITY

XII. M+W DATA

I. UNIT DESCRIPTION AND FUNCTION


Dental LED curing light M+W Superlite Power Pen/ Color Pen is a medical device - specialized light source for intraoral polymerization of dental materials, sensitive to the blue part of the light spectrum.

The device is intended for use only by a qualified dental practitioner and in a dental practice.

M+W Superlite Power Pen/ Colour Pen consists of a cordless handpiece and charging adapter. The Superlite Power Pen also includes a holder.

M+W Superlite Power Pen/ Color Pen is manufactured in conformity with the requirements of Regulation on medical devices MDR 2017/745 and standards ISO 13485:2016, ISO 9001:2015.

II. SYMBOLS

	Caution! Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Consult instructions for use Indicates the need for the user to consult the instructions for use.
	Dangerous voltage To indicate hazards arising dangerous voltage.
	Hazardous light emission To indicate hazards arising from light radiation.
	Hazardous thermal effects To indicate hazards arising from thermal effects.
	Manufacturer Indicates the medical device manufacturer.
	Date of manufacture Indicates the date when the medical device was manufactured
	Medical Device Indicates the item is a medical device
	Unique Device Identifier Indicates a carrier that contains Unique Device Identifier information
	Catalogue number Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Batch code Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Applied part type B according to electric safety classification
	Temperature limit Indicates the temperature limits to which the medical device can be safely exposed.
	Humidity limitation Indicates the range of humidity to which the medical device can be safely exposed.
	Waste Electrical and Electronic Equipment (WEEE) According to Directive 2012/19/EU, this symbol indicates that the product should not be disposed of as urban waste at the end of its operating life.
	Fragile Indicates a medical device that can be broken or damaged if not handled carefully.
	European Conformity Indicates conformity with local laws and regulations within the European Economic Area

III. SAFETY PRECAUTIONS

GENERAL WARNINGS

M+W Superlite Power Pen/ Color Pen is a Class I medical device and it meets the strict requirements of the Medical Devices Regulation - MDR (EU) 2017/745. In order to be used safely for staff and patients, the following rules must be observed:

- Do not allow unauthorized and untrained personnel to use the device to avoid risks.
- Disconnect the device from the mains after completing the procedures.
- Do not use or store the device in a dusty environment.
- Do not expose the device to direct sunlight.
- Do not spray disinfectant directly into the device - only rubbing with a swab drained of disinfectant is acceptable.
- Do not get wet or drop liquid on the device, cables, adapter to avoid electric shock or damage to the device.
- Store the device in a dry place, moisture can cause electric shock and damage.
- In case of a problem, disconnect the device from the mains, do not make attempts to repair, take the device to a service center.
- The device must not be used if any of its parameters are not normal (timer, light intensity, heat radiation).
- Strong electromagnetic fields in the building can cause interference and malfunction of the device. If their source cannot be determined, change the location of the device and plug it into another socket or other room, even in another building.
- Opening and repairing the appliance may only be carried out by authorized service technicians from the manufacturer.
- Only original M+W Superlite Power Pen/ Color Pen parts must be used when replacing defective parts. The warranty of the device does not cover the damage caused by the use of non-original spare parts. The device or any of its parts must not be disassembled while it is connected to the mains!

- Fragile! Use caution when transporting, using and storing the device!

- According to Directive 2012/19/EEC, this symbol indicates that the product should not be disposed as a general waste at the end of its lifespan. The product must be taken to a specialized center for the separate collection of electrical and electronic equipment according to local regulations. Proper disposal of equipment that is no longer used prevents negative consequences for the environment and human health!

- In accordance with the requirements of MDR (EU) 2017/745, user and / or the patient must report any serious accident that have occurred during use of the device to the manufacturer and the competent authority of the Member State in which the user/patient is established.
- All packaging materials of the product must be kept away from children to avoid risks of injury / suffocation.

SAFETY MEASURES AND RISKS

The device must be used in strict accordance with the Operating Instructions Manual.

1. Electrical safety

Before starting the appliance, make sure that the voltage and the type of plug correspond to the mains supply in the country. Use only the original adapter type SYS1638-0605-W2E.

Electrical safety is ensured by class II protection against electric shock according to EN 60601-1.

M+W Superlite Power Pen/ Color Pen must only be operated indoors, under the following conditions:

- temperature from + 10 ° to + 40 ° C;
- relative humidity 30 - 75%;
- lack of dust in the room;
- atmospheric pressure 700 - 1060 hPa;
- absence of chemically active and flammable substances;

- no part of the device should be wetted or immersed in water;
- the device or any of its parts must not be disassembled while it is connected to the mains!

Protect the cables of the appliance from insulation damage and breakage from sharp objects, strong pulling, rodents, chemicals. If such damage is noticed on the electrical cables, it is necessary to take the device immediately to the company service. The device must not be used with damaged cables.

In case of thunderstorms, the procedures must be stopped and the plug must be disconnected from the mains.

Risk: Failure to comply with these instructions may result in electric shock to users of the device.

2. Light radiation

M+W Superlite Power Pen/ Color Pen is a source of extremely intense light in the blue range, to which the human eye has a high sensitivity. This results in serious measures to be taken for patients, medical staff and accidentally nearby people, animals and plants.

As such, use protection goggles for the operator.

Irradiation of the eyes and skin with intense light carries a risk of damage from light and heat.

The light should never be directed at the eyes! Irradiation should be limited to the workplace area. The special safety goggles that meet the requirements must be used:

- to cover the eyes and temples tightly, even if the person is wearing optical glasses.
- be made of volumetric colored impact-resistant plastic.
- do not transmit light with a wavelength of 380 - 600 nm.
- reduce the intensity of the blue spectrum by more than 100 times.
- have a stable mechanical structure, no scratches, cracks and damage to its surface.

The device can be used only after a doctor's consultation on or by persons suffering from photo-biological reactions; persons taking photosensitive drugs; persons undergoing cataract surgery, persons with retinal diseases, etc.

The risk of improper irradiation is severe eye irritation, temporary spots in the visual field, severe visual impairment in direct radiation, to loss of vision.

3. Thermal radiation

The thermal effect is due to the absorption of the energy of the blue light in the tissues, during which the energy is converted into heat. The risk is only with prolonged overdose.

Risk of pain, burning of soft tissues.

4. Fire safety

- Keep the device away from solvents, flammable liquids and powerful heat sources.
- Do not expose to direct sunlight.
- Do not allow liquids and detergents to enter the device, as this may cause a short circuit and fire or cause potentially dangerous damage.
- If the product emits an odor or smoke – disconnect from the mains, do not attempt to repair it, take it to a service center.

Risk of fire, explosion and damage.

5. Contraindications

- The device can be used only after medical consultation on or by: persons with implanted cardiac pacemaker; persons suffering from photobiological reactions; persons taking photosensitive drugs; persons undergoing cataract surgery; persons with retinal diseases; people with allergies; people who have recently undergone cosmetic surgery on the face or lips, including injections of hyaluronic acid or botox; people with very sensitive skin or dermatitis, etc. If you are taking photosensitizers or medicines, check the package leaflet.


IV. TECHNICAL DATA

1. Operating voltage - 3.7V Li Ion battery type LIR 14500 (800mAh capacity)
2. Power consumption - mains adaptor: SYS1638-0605-W2E
 - 110-240VAC / 50-60Hz input voltage
 - Output: 5VDC / 1.2A micro USB socket
 - 0.5A max power consumption
3. Dimensions of cordless handpiece:
 - min. diameter – 18mm, max diameter - 25mm
 - length - 148 mm (without lightguide)
4. Weight of cordless handpiece (with lightguide and protection filter)
 - M+W Superlite Power Pen (plastic handpiece) - 80g
 - M+W Superlite Color Pen (metal handpiece) - 100 g
 - Holder – 276g
5. Polymerization mode – with soft start.
6. Light source – 2 band LED module with reflector optics.
7. Emitted light – blue, visible spectrum, 410-490 nm.
8. Light intensity - up to 1100 mW/sq.cm Intensity is constant and not dependent on battery charge level.
9. Emitting time: - 10 / 20 / 30 sec. / ±10%
10. Working ability after fully charged battery - up to 2100sec. (210 cycles x10sec).
11. Charge time - 3 hours for full charge and 2 hours to 80% charge.
12. Visual indication - color indication ring and audible signal every 10 sec.
13. Low battery warning. Fully discharged battery indication.
14. Fiberoptic lightguide (rotatable 360 degrees)
15. Working mode (1 min. work / 10 min. pause)
16. Degree of protection against electric current - applied part type B.

The manufacturer of this unit declares to provide on request all additional necessary technical documentation / information which will help the user's technical staff to service the parts of the unit which the manufacturer has claimed to be a subject for repair.

V. COMPLETE SET

1. Charge adaptor with micro USB output 5V DC / 1.2A - 1 pc.
2. Cordless handpiece - 1 pc.
3. Battery pack - 1 pc.
4. Lightguide - 1 pc.
5. Protection shield - 1 pc.
6. Protection cap - 1 pc.
7. Holder - 1pc.
(M+W Superlite Power Pen)
8. Operations guide - 1 pc.

Control and indication:

1. USB charging connector
2. Start / Stop button
3. Illuminated indicator ring
4. Battery pack

VI. PREPARATION AND SEQUENCE OF OPERATION

M+W Superlite Power Pen/ Color Pen is designed for use as a battery powered unit.


1. Take the curing light out of the packing and assemble the holder by inserting the upper part into the base part (as shown on picture below). Place the holder on a horizontal surface. The holder is only included with the Superlite Power Pen. For the Superlite Color Pen, it is optionally available in the colours black, white and silver.

2. Mount lightguide and protection shield or protection cap.

See VII. ROUTINE CARE AND MAINTENANCE for instructions of initial cleaning and mounting.

NOTE: Lightguide and protection cap must be sterilized before every use! They are delivered in a non-sterile state and must be sterilized before and first use. Before each patient, the light guide must be cleaned mechanically and then with a disinfectant and sterilised in an autoclave at 134 degrees. The temperature must be maintained for at least 5 minutes.

For polymerization of posterior teeth
(attach as shown)
(1) Protection cap

For polymerization of anterior teeth / general application

Connect the power adapter to charge battery if needed.

The device can be operated during charging.

Unplug the charger after charging the device or at the end of the day.

3. Indication

3.1. Color indication ring

Deactivated mode - Color indication ring is not lit.

Work mode – polymerization.

Battery charge / sound signal.

Fully charged battery – device can be unplugged from adapter.

Overheating – **Thermal protection activated**, wait until device cools down.

Battery is completely discharged - plug into the charging adapter.

3.2 LED-Modul

LED module indication:

Blinking means that the battery is discharged - there is only 20% charge left, you have the option of 5-25 starts of 10 seconds to full discharge.

Press the start and stop button once: 10 seconds of light polymerisation

Press the start and stop button twice: 20 seconds of light polymerisation

Press the start and stop button 3 times: 30 seconds of light polymerisation


After the set time has elapsed, the unit will turn off the output light.

If desired, the dentist can interrupt the cycle by pressing the Start / Stop button again.

During operation, the fiberoptic lightguide must be at a distance of at least 2 mm from the obturation

Diagram of curing mode

5. Overheating protection

The unit is equipped with overheating protection, which is activated if the temperature of the LED module rises up to 50°C. When overheating protection is activated the color ring will blink in yellow next 60sec, then is deactivated until the temperature of the handpiece reaches 35°C. It is necessary to wait for device's cooling down.

6. Battery – maintenance, charging and replacement

The device uses Li-Ion battery type LIR 14500 / 3,7V / 800 mAh.

Charging the battery - by connecting the charging adapter to the mains.

Attention! Connect properly to the socket!

Unplug the adapter after charging.

Charge adaptor is SYS1638-0605-W2E - 100-240VAC input and 5VDC output - micro USB type.

It is not allowed to use another type power supply.

Low Battery indication

- rare flashes of the output light from the lightguide - when 15-20% of the charge remains of the battery. It can be used for polymerization as the intensity is stabilized.

- With a fully discharged battery, the color ring will flash blue only a few times.

When charging - the light ring flashes slowly in blue, and when the battery is 100% charged - it lights up continuously.

Even if it remains switched on for charging for a longer period of time, there is no danger of overcharging or battery damage. After a full charge, the battery provides up to 2100 seconds work.

It is possible to work with the device during charging.

Monitor the battery charge and recharge if a low level is indicated.

Unplug the charger after charging the device or at the end of the day.

When to replace the battery

- at reduced capacity - frequent charging is required;

- in case of battery failure.

To replace a Li-Ion battery

- Disconnect the adapter from the socket and the USB socket.

- Place the handpiece on a non-slip surface and unscrew Part 3. Release Part 2.

Carefully remove the connector (4).



- The installation of a new battery is done in the reverse order.

- When assembling the handpiece - place Part 2 (containing the battery) on a non-slip surface and press Part 1 to Part 2 at the same time. At the same time, screw the cylindrical Part 3 back to the front conical Part 1.

Note: The parts of the assembly between 1 and 2 must fit!

Do not twist the battery cables!

Use only original Li-Ion batteries by the device manufacturer (LIR 14500 / 3,7V / 800 mAh)!

The battery must not be short circuited, mechanically damaged, immersed in water or to be heated or thrown into a fire.

- Non-functioning batteries can be taken to battery collection points or returned to the device manufacturer.

- When disposing of the battery, observe the environmental protection regulations!

VII. ROUTINE CARE AND MAINTENANCE

1. Examining of lightguide / protection cap.

The light guide is fitted by pushing it on and removed by pulling it off (can be rotated 360 degrees).

Every day the lightguide tip should be examined for:

- stuck restorative material (it should be removed very carefully with blunt non-metal object);
- mechanical damages on the lightguide.

Remove the light guide and looking through it at strong light: if there are darkened areas over 10% of the area, replace it with a new one.

Lightguide should be sterilized by autoclave at 134° C for 5 min.

The same must be done for protection cap.


The protective filter is wiped only with disinfectant.


2. Cleaning of the unit.

For disinfection of the unit, spray the disinfection agent onto a piece of soft cloth / cotton and clean the handpiece and power adapter.

Do not use abrasives or solvents as these may damage the plastic parts of the unit!

Do not spray directly onto handpiece or in button / USB charge connector!

VIII. WARRANTY CONDITIONS

1. Warranty period of M+W Superlite Power Pen/ Color Pen dental curing light is 24 (twenty four) months from the date of purchase. If date of purchase is not written, warranty starts at date of production.
2. During the warranty period the replacement of the defective parts is done free of charge by manufacturer's service. Lightguide is not covered by warranty. LED module and battery have 6 months warranty.
3. If during the guarantee period the unit is damaged, due to incorrect operation (mechanical, chemical, thermal or electrical damages), damage caused by improper use, storage or any other reasons in user's fault, the repair is paid by the user.
If the charging cable is damaged, it must be replaced. The device must not be charged or operated with a damaged cable.
- 3.1 If liquids, aggressive and flammable substances and their vapors have entered the device, the battery must be removed immediately and sent to the manufacturer's service where the damages are established. In case of such damage, the warranty is lost.
- 3.2 No damages or complaints are accepted as a result of electric shocks, thunderstorms, non-compliance with electrical safety measures or insufficient protection of patients, personnel and other persons from light radiation.
- 3.3 Claims due to improper or insufficient security and care, protection and security during transport, unpacking, relocation, operation and storage of the device are not accepted. The warranty is void for the above events.
- 3.4 The device must only be transported to manufacturer's service only in original packaging in order to avoid unwanted damage. The light guide must be packed separately from the cordless handpiece. The warranty card must be presented (this Guide, with completed Chapter XII M+W data) or a copy of the sales invoice showing the serial number of the device.
- 3.5 When carrying out repairs by unauthorized persons outside the company service or if non-original elements / parts are used, the user loses the right to free service.
4. The manufacturer does not owe compensation for lost profits in case of damage or incompleteness operation of the device, whatever the cause.
5. In case of disputes in connection with the application and interpretation of this Operating instruction manual, they will be decided by the court in Germany, by virtue of the current German legislation.
6. The warranty service is performed in the company service at the address:

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--

X. FAQ

Problem	Solution
No response by button starting	Fully discharged battery. Need to connect to charger to begin charge. Use only the original adapter type SYS1638-0605-W2E. If button continues not to respond - send the device to service.
The button does not respond properly	When working without barrier sleeves, the problem with the button may be the entry of composite, adhesive or disinfectant liquid when the device is used without protective sleeves. You need to clean with a cotton swab soaked in alcohol. If the button still does not respond correctly - send the device to service
Meaning of INDICATION ring colors:	<u>Flashing blue slowly</u> : in charging mode <u>constant blue</u> : battery 100% charged <u>Flashing yellow</u> : the lamp is overheated <u>infrequent flashing of the output light from the light guide</u> : the battery charge is at 15-20% <u>Flashing blue briefly</u> : the battery is discharged
Weak curing effect	It is necessary to clean the light exit window of the light guide with a suitable non-metallic object and to remove any adhering composite (restorative material) or other impurities with alcohol. If the result is poor - send the device to the service.
The patient feels discomfort during curing process especially during long curing time 20-30sec	Reduce the exposure time to 10 seconds and take a short break before the next exposure. The sensitivity is increased, especially for exposures in the near-pulse range.
Power cord or mains adaptor is damaged	Replace charging cable.
Entering some liquid into device	Send to service.
After drop on the floor device is not working properly	Send to service.

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 <p>TECHNICAL FILE LED curing light</p>	<p>Declaration of conformity</p> <p><i>Developed in conformity with MDR (EU) 2017/745</i></p>	TD 7.2
		<i>Revision 02</i>
		<i>Page 1 of 2</i>

Manufacturer: **BG LIGHT LTD**

SRN: BG-MF-000019812

Address: 155, Vasil Aprilov blvd., 4027 Plovdiv, Bulgaria
 Tel.: +359 32 644089, +359 888 809256
 BULSTAT 115841960, VAT N: BG115841960

DENTAL CURING LIGHTS

Basic UDI-DI: 3800501374200000VX

Intended use: M+W Superlite and A.M.Edelingsh (ORTHO) are LED curing lights designed to harden photocomposite materials used in dental practice.

N	Product code BGL	Product code M+W	Name of device	UDI-DI (GTIN)
1	200-003.01MW	200027	M+W Superlite Power Pen cordless (black)	3800501374405
2	200-003.02MW	200025	M+W Superlite Power Pen cordless (silver)	3800501374412
3	200-003.03MW	200026	M+W Superlite Power Pen cordless (white)	3800501374429
4	200-003p.04MW	200000	M+W Superlite Color Pen (white)	3800501374436
5	200-003p.05MW	200001	M+W Superlite Color Pen (orange)	3800501374443
6	200-003p.06MW	200002	M+W Superlite Color Pen (light green)	3800501374450
7	200-003p.09MW	200003	M+W Superlite Color Pen (blue)	3800501374467
8	200-003p.10MW	200004	M+W Superlite Color Pen (violet)	3800501374474
9	200-003p.12MW	200005	M+W Superlite Color Pen (pink)	3800501374481
10	200-003p.02MW	200006	M+W Superlite Color Pen (dark grey)	3800501374498
11	200-004MW	100045	M+W Superlite 1300	3800501374719
12	200-004h MW		Handpiece for M+W Superlite 1300	3800501374856
13	200-009.01 MW	144000	A.M.Edelingsh (black)	3800501374658
14	200-009.03 MW	144006	A.M.Edelingsh (white)	3800501374665
15	200-009.01ort MW	144010	A.M.Edelingsh ORTHO (black)	3800501374832
16	200-009.03ort MW	144011	A.M.Edelingsh ORTHO (white)	3800501374849

SPARE PARTS FOR DENTAL CURING LIGHTS

Basic UDI-DI: 3800235341000003G

Intended use: Spare parts for M+W Superlite and A.M.Edelingsh - LED curing lights designed to harden photocomposite materials used in dental practice.

N	Product code BGL	Product code M+W	Name of spare part	UDI-DI (GTIN)
17	200-003H.01MW	200030	Holder for M+W Superlite Power Pen (black)	3800235340325
18	200-003H.02MW	200028	Holder for M+W Superlite Power Pen (silver)	3800235340332
19	200-003H.03MW	200029	Holder for M+W Superlite Power Pen (white)	3800235340349
20	200-003b.01MW	100036	Battery pack for M+W Superlite Power Pen (black)	3800235340356
21	200-003b.02MW	100091	Battery pack for M+W Superlite Power Pen (silver)	3800235340363
22	200-003b.03MW	100035	Battery pack for M+W Superlite Power Pen (white)	3800235340370
23	200-003pb.04MW	200007	Battery pack for M+W Superlite Color Pen (white)	3800235340387
24	200-003pb.05MW	200008	Battery pack for M+W Superlite Color Pen (orange)	3800235340394
25	200-003pb.06MW	200009	Battery pack for M+W Superlite Color Pen (light green)	3800235340400
26	200-003pb.09MW	200010	Battery pack for M+W Superlite Color Pen (dark blue)	3800235340417
27	200-003pb.10MW	200011	Battery pack for M+W Superlite Color Pen (violet)	3800235340424
28	200-003pb.12MW	200012	Battery pack for M+W Superlite Color Pen (pink)	3800235340431
29	200-003pb.02MW	200013	Battery pack for M+W Superlite Color Pen (dark grey)	3800235340448
30	200-004.001MW	100047	Battery for M+W Superlite 1300	3800235340455

31	200-009.003.01MW	144002	Battery pack for A.M.Edelingsh (black)	3800235340462
32	200-009.003.03MW	144004	Battery pack for A.M.Edelingsh (white)	3800235340479
33	200-003.004MW	100092	Power adapter for M+W Superlite Power Pen / Color Pen	3800235340486
34	200-004.004MW	100049	Power adapter for M+W Superlite 1300	3800235340493
35	200-009.004MW	144003	Power adapter for A.M.Edelingsh	3800235340509
36	200-009-001MW	144001	Barrier sleeves for A.M.Edelingsh	3800235340516
37	102-001MW	135408	Eye protection shield (red) for M+W Superlite	3800235340547
38	100-009MW	100046	Bleaching lightguide for M+W Superlite 1300	3800235340554
39	100-001MW	100048	Fiberoptic lightguide for M+W Superlite	3800235340561
40	100-001-1MW	100043	Special lightguide (light intensity -30%) for M+W Superlite	3800235340578
41	102-002MW	100094	Protection silicone cap (grey) for M+W Superlite	3800235340585

Classification: Active device of Class I of the Regulation on medical devices - MDR (EU) 2017/745

Classification is done by the manufacturer according to Regulation on medical devices - MDR (EU) 2017/745, Rule 13, Annex VIII.

Notified body: TUV NORD Polska Sp. z o.o., ul.Mickiewicza 29, 40-085 Katowice, Poland. NB number: 2274.

This Declaration of conformity is valid only in combination with our certificates of Notified body TUV NORD Polska Sp. z o.o. Certificates N: AC090 100/1971/4047/2020, AC090 MD/1971/4047/2020.

The manufacturer declares under its sole responsibility that the products are developed and produced in conformity with MDR (EU) 2017/745 and the following applicable standards:

EN 60601-1:2006 +AC:2010+A1:2013+A12:2014	Medical electrical equipment - Part 1: General requirements for safety
EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests
EN 62304:2018	Medical device software. Software life cycle processes.
EN 62353:2015	Medical electrical equipment. Recurrent test and test after repair of medical electrical equipment
EN 62366-1:2015+AC:2016	Medical devices. Application of usability engineering to medical devices.
EN ISO 14155:2011+AC:2011	Clinical investigation of medical devices for human subjects - Good clinical practice
EN ISO 14971:2019Q	Medical devices – Application of risk management to medical devices.
CEN ISO/TR 24971:2020	Medical devices - Guidance on the application of ISO 14971 (ISO/TR 24971:2020)
EN ISO 15223-1:2016	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
EN ISO 20417:2021	Medical devices — Information to be supplied by the manufacturer
EN ISO 10993-1:2018	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
Directive 2012/19/EC	Directive on waste electrical and electronic equipment (WEEE)

All company products are manufactured under the current Quality Management System, ISO 9001:2015 and ISO 13485:2016.

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M+W Dental  Müller & Weygandt GmbH · Reichardsweide 40
BG LIGHT LTD · 155, Vasil Aprilov blvd. · BGR - 4027 Plovdiv · tel.: +359 32 644089