// Pressure Moulding Technique

BIOPLAST® XTREME PRO Kit

Herstellanleitung / Instruction manual /
Mode d'emploi / Manual de instrucciones

SCHEU
Dental Technology

Materialliste

BIOPLAST® XTREME 4,0 mm

BIOPLAST® 2,0 mm

DURAN®-Folie 0,75 mm

BIOPLAST®-Isoliermittel

Heiz-Support

Finier-Set

Granulatabdeckschablone

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BIOPLAST® XTREME 4,0 mm

BIOPLAST® 2,0 mm

DURAN® 0,75 mm

BIOPLAST® agent isolant

Support de chauffe

Kit de finition

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List of materials

BIOPLAST® XTREME 4.0 mm

BIOPLAST® 2.0 mm

DURAN® 0.75 mm

BIOPLAST® insulating agent

Heating Support

Finishing set

Pellet cover

Lista de materiales

BIOPLAST® XTREME 4,0 mm

BIOPLAST® 2,0 mm

DURAN® 0,75 mm

BIOPLAST® líquido separador

Soporte de calentamiento

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Cobertura de granulado

BIOPLAST® XTREME PRO Kit

Herstellanleitung	4 - 10
Instruction manual	11 - 17
Mode d'emploi	18 - 24
Manual de instrucciones	25 - 31
Heizzeiten und Codes	
Heating times and codes	
Temps de chauffe et codes	
Tiempos de calentamiento y códigos	32 - 33

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BIOPLAST® XTREME PRO

Instruction manual

Apply one thin layer of BIOPLAST® insulating agent only to those parts of the model that will be covered by the mouth guard. Pay attention that the insulation does neither dry nor accumulate in the recesses. Trim model flat and place it on model platform or into the pellets (using the pellet cover) with the clear V-shaped part of the BIOPLAST® XTREME material covering the desired extension of the anterior section (e.g. from 3-3). Make sure the side of the label with the imprint is facing the infrared heater when fixing the material. Heat the foil by setting the code or recommended heating time according to the instructions. Close pressure chamber and open once cooling phase has elapsed.

Instruction manual

GB

Cut material with scissors or scalpel on the model bottom and insulate the anterior section with CETRON® spray.

For the fabrication of a protection insert: embed the model vertically in the pellets and heat DURAN® 0.75 mm by setting the code or recommended heating time. Lock pressure chamber and open once cooling phase has elapsed.

Cut the BIOPLAST® XTREME to the desired extensions. Make sure the material is being cut straight edge, as cutting it diagonally might impair the optimum lamination with the 2nd BIOPLAST®. Cut the DURAN®, too, with scissors to the desired extensions.

Make sure the extension of DURAN® towards the alveolar ridge is some millimeters shorter than the extension of BIOPLAST® XTREME underneath, in order to avoid problems in laminating the material. The impression of the counter bite has to be considered, too.

Pierce holes into the BIOPLAST® XTREME – occlusally and interdentally – by means of a bur, in order to avoid air bubbles under the protection insert and between the BIOPLAST® material.

Clean BIOPLAST® XTREME with moist tissue paper to remove the insulation. Do not use any re-fatting cleaning agents.

Before pressurising the BIOPLAST® 2.0 mm, place the protection insert made of DURAN® in the front, the front of the model facing the pressure chamber.

Then heat BIOPLAST® 2.0 mm by setting the code or recommended heating time, lock the pressure chamber and open once cooling phase has elapsed. Cut material at model base.

Cover the margins with a damp cloth as heat protection. Heat the occlusal surfaces under the infrared heater using the heating support with the pressure chamber being closed, but not locked. Insulate occlusally the antagonist model with BIOPLAST® insulating agent, in order to avoid plaster residues. Then put models in the articulator for the impression of the counter bite. Cool down with compressed air and water and remove the models.

Cut the material with scissors to desired extensions.

Polish the edges with BIOPERM® trimmer or DIMO®PRO wheel. Cut the mouth guard diagonally to remove sharp edges in the labial and palatal area. Finish the rims using the OSAMU polisher at a low RPM exerting low pressure.

Die tabellarisch aufgeführten Ergebnisse gelten für alle Farbvarianten der Folientypen BIOPLAST® XTREME und BIOPLAST® XTREME PRO.
The heating times and codes given in the table above are valid for all colour variations of BIOPLAST® XTREME and BIOPLAST® XTREME PRO.
Les résultats énumérés dans ce tableau sont valables pour toutes les couleurs de BIOPLAST® XTREME et BIOPLAST® XTREME PRO.
Todos los tiempos de calentamiento indicados en esta tabla son válidos para todos los colores de BIOPLAST® XTREME y BIOPLAST® XTREME PRO.

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**Three Mouthguard sports
protection types –**
meeting all requirements.

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All mouthguard types as described above consist at least of 2 layers, guaranteeing maximum safety and comfort.

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SCHEU-DENTAL GmbH