# FREEPRINT® Ortho

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# FREEPRINT® Ortho

Light-curing resin based on (meth)acrylate, biocompatible, for the generative fabrication of orthodontic bases, occlusal orthos, drilling templates, X-ray templates

for DLP printers with LED with 405 nm / UV-LED 378 - 388 nm

## Important notes

This is a medical device, only to be used by trained specialist personnel.

## Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 2).
- ▶ Minimum material thickness for the design is 1.5 mm
- ▶ Maximum curing depth\* at direct post-exposure:  
clear: 3,5 mm  
\*In case of large objects and exposure on both sides, the material thickness can be up to 7 mm (Example FREEPRINT® ortho – with a curing depth of 3,5 mm).  
Do not use heat-based methods for disinfection or sterilisation. This could possibly deform the workpiece.
- ▶ Polish surface mechanically.
- ▶ Processing temperature 23°C ± 2°C.

## Safety

- ▶ Please follow the instructions on the safety data sheet!
- ▶ Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- ▶ Avoid direct contact with the liquid material and the components prior to post-curing. Irritating to eyes and skin (sensitisation is possible).
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.

## Storage

- ▶ **FREEPRINT® ortho** is to be stored dry (at 15°C - 28°C) and protected from light. Minimal influence of light can already induce polymerisation.

## General

- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.

## Contraindication

Contains (meth)acrylics

Some ingredients of **FREEPRINT® ortho** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® ortho** only insert intraorally in completely polymerised state.

## Adverse effects

Product may cause allergic reactions.

**FREEPRINT® ortho** was developed for use in dental laboratory and must be used in accordance with the instructions for processing and safety. DETAX will not be responsible for damages caused by faulty or improper use of system and materials.

## Indication:

Orthodontic bases, drilling and X-ray templates for implant dentistry, occlusal orthos, fixation and transfer keys

## Processing:

at 23°C ± 2°C

## Storage

## Ordering information:

**FREEPRINT® ortho**  
1.000 g bottle, clear **03988**

**FREEPRINT® ortho UV**  
1.000 g bottle, clear **03989**

**FREEPRINT® splint**  
1.000 g bottle, clear **03066**

**FREEPRINT® splint UV**  
1.000 g bottle, clear **03709**

**FREEPRINT® model**  
1.000 g bottle  
sand **03065**  
ivory **03779**  
grey **03781**

**FREEPRINT® model UV**  
1.000 g bottle  
sand **03778**  
ivory **03780**  
grey **03782**

**FREEPRINT® cast D**  
500 g bottle, red **02890**

**FREEPRINT® cast UV D**  
500 g bottle, red **03710**

More information at  
[www.detax.de](http://www.detax.de)

# DETAX

GmbH & Co. KG  
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany  
Telefon: 07243/510-0 Fax: 07243/510-100  
[www.detax.de](http://www.detax.de) · [post@detax.de](mailto:post@detax.de)

Made in  
Germany

06/2016

# FREEPRINT® Ortho

## Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

## Construction process

Generation of a Print Job complying with machine and material parameters

## Post-processing

If possible, post-processing should commence immediately following the construction process. After raising the platform, a drip time of approx. 10 minutes is recommended.

## Pre-cleaning

Remove construction components from the platform and clean in a separate vessel with isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath.

## Cleaning

Then thoroughly clean the openings, cavities and gap areas, if necessary also with compressed air, and, if applicable, remove the construction components carefully from the support structure.

## Main cleaning process

The main cleaning process is performed in a separate vessel with fresh isopropyl alcohol (purity  $\geq 98\%$ ) for 3 min. in an ultrasonic bath. Prior to drying, check the openings and additional bore holes for residues.

## Drying

Heat the construction components for 30 min. in an oven to approx. 40°C to remove the solvent residues from the cleaning process.

## Post-exposure

Post-exposure is performed with a xenon photoflash unit Otoflash G171 with 2 x 2000 flashes under inert gas conditions (nitrogen), rotate components in between.

## Surface processing

Polish surface mechanically