IP LIGHTING PRO

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■ 1 General information

The illustrations, descriptions and properties in this booklet are not binding. Apart from the main characteristics, the manufacturer reserves the right to make changes at any time due to, to make changes at any time due to technical or commercial requirements. The device has been designed and realised for use in the field of dental technology and dentistry. Other applications are not intended by the manufacturer. He consequently declines all responsibility for any damage resulting from improper use of the device. Improper use of the device. Modifications, unauthorised interventions (even of a minor nature) and the use of non-original spare parts, original spare parts release the manufacturer from any responsibility and forfeit the right to warranty. The appliance may only be operated by adult and appropriately trained personnel.

■ 2 Electrical connection

Before connecting the appliance to the power supply, please ensure that the data on the rating plate matches that of your power supply and that an earthing cable is available. Insert the plug of the loose cable coming from the appliance into an appropriate socket. The main switch of the appliance must be switched off. The power supply to the appliance must comply with the instructions from IP Division GmbH. However, the latter is not responsible for the electrical connection on site. The supplied cable is of type H05 with a diameter corresponding to the maximum current consumption (see technical data). The power cable has a length of 2 metres. When using extension cables, please note that

- the average of the cable corresponds to the length. The average must correspond approximately to that of the cable supplied.
- the extension cable does not present a tripping hazard.
- the connection to the electrical power supply complies with all local safety regulations.

Warning of dangerous electrical voltage! Warning signs are safety signs warning signs that warn of risk or danger. DIN 4844-2 Safety labelling Part 2: Presentation of safety signs. BGV A 8: Safety and health signs at the workplace safety and health signs at the workplace - Appendix 2 W08 Warning of dangerous electrical electrical voltage (lightning arrow warning sign).

3 Device structure & scope of delivery

- 1 Main switch
- 2 Timer
- 3 Door Handle Polymerisation chamber
- 4 Turntable
- 5 Minutes key
- 6 Seconds key
- 7 Start / Stop
- 8 IEC socket with microfuse
- 9 Power cable

4 Safety regulations

WARNING: To avoid the risk of electric shock, injury and fire when using electrical tools, the following instructions must be read carefully and observed at all times:

- A clearly organised workplace helps to avoid accidents.
- Environmental influences: Electrical appliances must never be exposed to rain or operated in wet and damp environments. Ensure adequate lighting. Never use the appliance in the vicinity of flammable liquids or gas.
- Avoid direct contact with earthed objects, e.g. tubes, spotlights, refrigerators.
- Keep children away! Third parties are prohibited from handling the appliance and touching the cable, as well as from standing in the immediate vicinity of the appliance.
- Set up the appliance in a safe, dry place that is inaccessible to children.
- Do not overload the appliance! Observe the power ratings specified on the rating plate.
- Always wear UV safety goggles.
- Take care of the cable: Never pull on the cable, e.g. to lift the appliance or to pull the mains plug out of the socket. Protect the cable from high temperatures, oil and sharp edges.
- Avoid adopting an unsafe working posture.
- Look after your instrument carefully. To ensure safe operation, keep your appliance clean, check the cables regularly and have them replaced by qualified specialists if in doubt. The same applies to extension cables. The buttons and handles of the appliance must be clean and dry.
- Ensure that there are no loose objects such as tool spanners on the appliance when it is being operated.
- To avoid accidental operation of the appliance, make sure that the appliance is switched off when you plug in the mains plug.
- The appliance may only be operated in closed and well-lit rooms.
- Danger due to distraction: Lack of attention can lead to accidents. We recommend that third parties keep their distance from the appliance and its operator.
- Check for damage: Before using the appliance, moving parts and safety devices must be thoroughly checked for functional safety. Make sure that no parts are jammed or damaged, that all components are correctly installed and that all factors affecting the work process are optimised. Damaged parts must be replaced or repaired by technical customer service, unless other instructions are expressly given in the operating instructions.
- Caution! For your own safety, only use tools or accessories that are listed in the operating instructions or suitable technical literature.
- Spare parts may only be replaced by qualified personnel. This electrically operated appliance complies with the applicable safety regulations. Repairs may only be carried out by appropriately trained customer service personnel.

5 Start of operation

Your device is packed in a cardboard box. Carefully remove the adhesive tape and open the box. Lift the device out with both hands and place it on a stable work surface. Remove the polystyrene packaging and the plastic film.

5.1 Insert

- Open the door (3) and place the model in the centre of the turntable.
- Do not insert any objects that could hit the walls of the chamber during the rotation.
- Carefully close the door with the safety microswitch.

5.2 Switching on

After closing the door (3) and pressing the main switch (1), the duration of the operating cycle can be set using the two buttons minutes (5) and seconds (6). The operating temperature is approx. 45-50°C.

5.3 Replacing the lamps

The device contains four different LED modules with a total of 48 high-performance LEDs that cover the full light spectrum from 380 to 650 nm. The LED lamps should never be replaced by the user. They have a very long service life. However, should a defect occur, please contact an authorised customer service centre.

5.4 Tested plastics and polymerisation times

The following veneering composites have been successfully tested in cooperation with master laboratories (as of April 2022, the list is constantly updated):

Manufacturers	Series
Ivoclar	Nexco
Heraeus-Kulzer	Signum
GC	Gradia
3M	Sinfony
Shofu	Ceramage
Schütz-Dental	Dialog
Bredent	Crea.Lign

Recommended polymerisation times

We recommend the following polymerisation times for the initial application:

Plastic type	Impression tray	Composite	3D-Resin
Time	40-90 seconds (depending on layer thickness)	2-3 minutes (depending on layer thickness)	3-5 minutes (depending on layer thickness)

These are general guide values for first-time users. The times may vary depending on the layer thickness, manufacturer and desired hardness.

6 Care

Use a slightly damp cloth to clean the appliance housing. Ensure that no moisture penetrates the openings and do not use any solvents or abrasive cleaners.

6.1 Place of use

The room should be well ventilated, dry and clean.
Indoor temperature: 0°C - 35°C.
Max. relative humidity: 85%.
Switch off the appliance immediately after use.

6.2 Warning notes

Do not spill any liquids on the appliance.
Observe the recommended polymerisation cycles - longer polymerisation times can have a negative effect on the final hardness of the materials.

7 Technical data

Dimensions	28 x 24 x 29 cm (L x W x H)
Turntable dimensions	160 mm
Mains voltage	220V, 50/60 Hz
Power	115 W
Wave length	380-650 nm
Timer	Electronic timer from 0-99 min
Weight	10 kg

8 Troubleshooting

All of the following measures may only be carried out when the appliance is disconnected from the power supply. Do not carry out any maintenance work, adjustments or repairs that are not mentioned in these operating instructions. If necessary, please call our technical customer service.

Defect	Cause	Remedy
The device does not start.	Fine-wire fuse in the device defective	Check the fuse and replace it if necessary
	Cable defective or loose	Check that the respective plugs of the connection cable are correctly seated on the appliance and on the power source. Insert a new cable.
	Main switch is defective	Contact customer service
The LED lamps do not light up.	Power supply to the lamps is interrupted	Check that the respective plugs of the connection cable are correctly seated on the device and on the power source.
The LEDs do not switch on despite the cable being connected correctly	Door is not closed	Close the door completely
	Door safety micro switch defective	Contact customer service
	Light source defective	
The fan does not work during operation	Door is not closed	Close the door completely
	Defective fan	Contact customer service
Turntable does not rotate	Door is not closed	Close the door completely
	Turntable defective	Contact customer service

9 EC - Declaration of Conformity

The following appliance has been tested in accordance with Directive 89392/EEC and subsequent versions. Conformity is confirmed by the type plate with the CE marking on the appliance and by the declaration of conformity enclosed in the appendix.

The manufacturer
IP Dental Division GmbH
Dachauer Str. 95a
85778 Haimhausen
Deutschland

hereby declares that the following device fulfils all the requirements of the following directives and standards:

Produktbezeichnung:	IP Lighting Pro
Art. Nr.:	6340
2006/42/EC	Machinery Directive
2014/30/EC	Electromagnetic compatibility (EMC)
2014/35/EC	Low Voltage Directive
2012/19/EC	Disposal of waste electrical & electronic equipment
EN 61000-6-3	Interference emission
EN 61000-6-1	Interference resistance
EN 60335-1	Machinery safety
EN 61326-1	Electrical measurement, control and laboratory equipment
EN 61010-1	Safety regulations for electrical measuring, control and laboratory equipment

Haimhausen, 20.03.2024

IP Dental Division GmbH

10 Disposal in accordance with EU directives

- This appliance must not be disposed of with household waste!
- Disposal exclusively by certified specialised companies
- Registration number Germany 21536414

■ 11 Warranty and disclaimer

P Division GmbH declares that the IP Lighting Pro product is free from defects in material and workmanship.

Warranty after date of purchase:

- **12 months on the entire device.**

The warranty is only valid if the device was purchased undamaged from the manufacturer IP Division GmbH or authorised retailers.

It does not apply in the event of improper use of the device or if modifications, unauthorised interventions or repairs have been carried out by persons who are not authorised service technicians of IP Division GmbH.

If the product is found to be defective in material or workmanship within one year of the date of purchase, it must be returned to IP Division GmbH or an authorised service centre at the sender's expense with the name, address and proof of date of purchase, together with a brief description of the defect.

IP Division GmbH will then, at its own discretion, repair or replace the defective device at no cost to the purchaser. The warranty for possible further repairs or measures remains valid for the remainder of the warranty period as described above. The warranty only covers the repair and maintenance of the defective product. There is no further express or implied warranty. IP Division GmbH is not responsible for incidental or consequential damages or the resulting cost of repair or replacement as described above.

IP Dental Division GmbH
Dachauer Str. 95a
85778 Haimhausen
Germany
Tel.: +49 - 81 33 - 91 79 480
E-Mail: info@ipdent.de
Web: www.ipdent.de