Permaplast LH Hybrid Spritze

Stand:07/1024

D Produktbeschreibung

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GB Product description

M+W Permaplast LH Hybrid Syringe is a light-curing filling material for restorations of all cavity classes.

Indications
Loss of tooth substance due to caries, trauma, tooth attrition/abrasion or development disorders.

Contraindications
• Do not use directly on the exposed pulp.
• Do not use the material where there are contact allergies or known allergies to any of the components.

Intended use
Filling material for restorations of all classes of cavities

Limitations of use
Do not use the material if a relatively dry working area or the recommended application technology is not possible.

Intended users
Dentist

Patient target group
People treated in the course of a dental procedure

Recommended use
1. Clean the tooth. Remove any cleaning residue with a water spray. Determine the VITA shade.

Note: The use of a rubber dam is recommended.

- Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Protect areas close to the pulp.
- If necessary, etch using your preferred etching technique (e.g., Total-Etch). Observe the manufacturer's instructions for the material used.
- Apply a bonding agent in accordance with the manufacturer's instructions.
- If necessary, apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.
- Apply the desired shade of the M+W Permaplast LH Hybrid Syringe with a maximum layer thickness of 2 mm into the cavity. Light-cure individual layers for 40 seconds each at a device light intensity of $\geq 600 \text{ mW/cm}^2$ or for 20 seconds at a device light intensity of $\geq 1000 \text{ mW/cm}^2$.
- Removal of excess material, contouring and fine processing are done, for example, with diamond finishing tools and flexible polishing discs of progressively finer grit sizes. M+W Permaplast LH Hybrid Syringe can be given a high shine using polishing paste together with polishing cups or discs.

Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons. Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 600 mW/cm^2 . Place the light-curing unit as close as possible to the material.
- Do not remove the oxygen-inhibited layer if a further layer of filling material is to be applied, because this is required for the bond with the next layer.
- Follow the manufacturer's instructions for any other materials that are used with Universal hybrid composite.

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid materials containing eugenol, moisture and air containing oil!
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components in the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

The product complies with ISO 4049:2019, Type 1, Class 2, Group 1.
Polymerization depth $\geq 2 \text{ mm}$ with a 40 second light-curing time when the light intensity of the device is $\geq 600 \text{ mW/cm}^2$. Polymerization depth $\geq 2 \text{ mm}$ with a 20 second light-curing time when the light intensity of the device is $\geq 1000 \text{ mW/cm}^2$.

Storage/disposal

- Store at room temperature ($15 - 25^\circ \text{C} / 59 - 77^\circ \text{F}$).
- Close syringe immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date

Composition

Dental glass, polymer filling material, zeolite, SiO_2 , bis-GMA, UDMA, HDDMA, additives.
Inorganic filling materials: approx. 68 vol% (0.02 – 3.6 μm).

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

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