

Inhalt/Content:

20 kg



Stand:030423

SIMPLEE®

All you need. No more. But no less.

SIMPLE ROCK

Kunststoffveredelter Klasse IV Gips zur Herstellung von Zahnkranz- und Modellguss-Modellen. Thixotrope Konsistenz mit niedriger Expansion.

© Plastic-enhanced Class IV plaster for the fabrication of denture and model casting moulds. Thixotropic consistency with low expansion.

REF 700 910 goldbraun/golden brown
UFI: 3G4U-Q77K-ME9N-UWF6

REF 700 912 pastellgelb/pastel yellow
UFI: 1K4U-67WY-XE95-H818

Anmischverhältnis/Mixing ratio: 100 g: 20 ml dest. Wasser/
Distilled water

Anmischzeiten/Mixing times: manuell/Manually 60 Sek./sec.
Vakuüm/Vacuum 30 Sek./sec.

Verarbeitungszeit/Processing time: 5 - 6 min.

Abbindezeit/Setting time: 9 - 10:30 min.

Enformbarkeit/Demouldability: ca. 30 min.

Druckfestigkeit nach/Compression strength after 1h: 70 MPa

Expansion nach/after 2h: 0,08%

Expansion nach/after 24h: 0,09%

 Lifco Dental AB; Verkärestaregatan 1 · 745 39 Enköping, Sweden