

Inhalt/Content:

20 kg



Stand:030423

SIMPLEE®

All you need. No more. But no less.

SIMPLE FLOW

Dünflüssiger Klasse IV Sockelgips. Gute Fließfähigkeit bei niedriger Expansion und glatter Oberfläche. Splitcast und Modellssockel bleiben blasenfrei und passgenau.

© Low-viscosity Class IV base plaster. Good flowability with low expansion and smooth surface. Splitcast and model base remain bubble-free with a perfect fit.



REF 700 913 extraweiß/extra white

UFI: 8XWD-PC9V-OE9G-C281



REF 700 915 gelb/yellow

UFI: 1PWD-5C7P-3E90-D1GU

Anmischverhältnis/Mixing ratio: 100 g: 23 ml dest. Wasser/
Distilled water

Anmischzeiten/Mixing times: manuell/Manually 60 Sek./sec.
Vakuum/Vacuum 30 Sek./sec.

Verarbeitungszeit/Processing time: 5 min.

Abbindezeit/Setting time: 10:30 - 11:40 min.

Enformbarkeit/Demouldability: ca. 30 min.

Druckfestigkeit nach/Compression strength after 1h: 50 MPa

Expansion nach/after 2h: 0,06%

Expansion nach/after 24h: 0,06%

 Lifco Dental AB; Verkärestaregatan 1 · 745 39 Enköping, Sweden